

# Isoelectric Focusing

Craig Nelson

Manager Electrophoresis Technologies

BioDiagnostics Inc.

# Isoelectric Focusing

- Electrophoresis in a pH gradient
- Separation method that resolves protein markers on the basis of their isoelectric points
- Successfully used in clinical, research, and agricultural fields
  - Blood
  - Serum
  - Muscle Extracts
  - Seed Extracts

# Isoelectric Focusing

- Attributes of IEF
  - Ability to separate on total protein marker profiles or enzyme protein marker profiles
  - Proteins with same molecular weights will separate out by pH
  - Testing time is quick (~2- 3days)
  - Ability to have permanent record of the gel
  - Pre cast gels available
  - ┆ Versatile to a large number of crops for the agriculture industry
  - ┆ Seed or tissue extraction

# Isoelectric Focusing Sample Preparation

- Load seed or tissues into a 48 or 96-well microplate
- Utilize a seed crusher or cutter to extract the seed proteins




# Isoelectric Focusing Sample Preparation

- A specific quantity of extraction solution is added into the wells of the microplate with a multi-channel pipette
- Incubate samples overnight at 4C
- Centrifuge samples with the resulting solution ready to be applied to an IEF gel
- The microplate containing the sample solution can be frozen for sample retesting or future evaluation

# Isoelectric Focusing of Samples

- Prepare Multiphor unit
  - Turn on circulating water bath 12-15C
  - Level unit
  - Clean cooling plate




# Isoelectric Focusing of Samples

- Place 2 mL of water onto the center of the cooling plate
- Mark a reference point on the edge of gel
- Place the gel on the cooling plate, distribute the bead of water being careful to not trap any air bubbles
- Blot excess water from edge of gel
- Use blotting paper to remove moisture from the gel



# Isoelectric Focusing of Samples

- Prepare three precut electrode wicks
- Saturate two wicks with anode solution
- Place each anode wick parallel along the long edges of the gel





# Isoelectric Focusing of Samples

- Saturate the remaining wick with cathode solution
- Place the cathode wick along the center of the gel between the two anode wicks
- Prefocus the gel if this step is required
  - Place the electrode wires directly on the anode and cathode wicks
  - Lower electrode cover connect to power source, run electric current at specific wattage for 10 minutes

# Isoelectric Focusing of Samples

- Place paper or plastic template on each side of the gel at its designated position (1cm from the anode wicks)




# Isoelectric Focusing of Samples

- Use application tool to take off frame of paper template



# Isoelectric Focusing of Samples

- Final layout of the gel



# Isoelectric Focusing of Samples


- Pipette the prescribed quantity of sample solution into each sample well





# Isoelectric Focusing of Samples

- Focusing
  - Electrode wires are placed onto wicks
  - Safety cover is placed on unit
  - Connect electrode wires to power supply
  - Gel is run at constant voltage for a determined amount of time



# Isoelectric Focusing Silver Staining

- Place the recently run gel into a staining dish, pour 200 ml of fixative solution onto the gel, rock for 15 minutes
- Wash the gel while rocking for one hour in 500 ml of deionized water, repeat wash step with fresh water
- Allow the gel to completely dry overnight at room temperature



# Isoelectric Focusing Silver Staining

- Clean all staining glassware and stir bars with reducing wash solution
- Formulate the silver staining solution
- Prepare gel for staining by washing for 5 minutes in 200 ml deionized water
- Pour silver stain mixture over the gel's surface and begin rocking the gel
- Allow the gel to stain until bands reach desired intensity

# Isoelectric Focusing Silver Staining

- Discard the staining solution and rinse gel in deionized water
- Pour stop solution onto the gel, soak for 10 minutes
- Place the gel in 200 mL of deionized water, soak for 10 minutes
- Allow gel to air dry

# Isoelectric Focusing Enzyme Staining

- Thirty minutes before the completion of the gel run prepare the enzyme stain
- Pour the stain over the gel surface immediately after focusing
- Rock and heat the gel during staining




# Isoelectric Focusing Enzyme Staining


- Allow the staining to continue until the bands reach strong intensities
- Stop the reaction with 200 mL enzyme stop solution for 3 minutes
- Wash the gel with two cycles of 500 mL deionized water for one hour
- Allow the gel to dry
- The gel can now be analyzed and interpreted

# Analysis of Protein Markers

**Coomassie Blue Stain of Wheat**




**Silver Stain of Soybeans**




# Analysis of Protein Markers, cont.

**Silver Stain of Peppers**




**Alcohol Dehydrogenase of Tomatoes**






# Analysis of Protein Markers, cont.

**Silver Stain of Catfish Whiskers**



**Acid Phosphatase of Canola**





# Analysis of Protein Markers, cont.



# Applications of Isoelectric Focusing

- Agricultural Industry
  - Genetic Purity
  - Variety identification, including trait variety identification.
  - Parental Line Maintenance
  - Breeding Programs

Questions?